Rematrix Ortho

- GB** Acellular dermal matrix for soft tissue reconstruction
- DE** Azelluläre dermale Matrix für die Weichgewebe-Rekonstruktion
- FR** Matrice dermique acellulaire pour la reconstruction de tissus mous
- ES** Matriz dérmica acelular para la reconstrucción de tejido blando
- IT** Matrice dermica acellulare per la ricostruzione di tessuti molli

REF	Size/Thickness [mm]	Content per Box
RMT0505E	50 x 50 / 1,4 +/- 0,2	1 unit
RMT0510E	50 x 100/ 1,4 +/- 0,2	1 unit
RMT1020E	100 x 200/ 1,4 +/- 0,2	1 unit
RMT0505F	50 x 50/ 0,8 +/- 0,2	1 unit
RMT0510F	50 x 100/ 0,8 +/- 0,2	1 unit
RMT1020F	100 x 200/ 0,8 +/- 0,2	1 unit

 **MBP Medical Biomaterial Products GmbH**
Lederstraße 7
D-19306 Neustadt-Glewe
info@mbp-gmbh.de/ www.mbp-gmbh.de

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Rematrix Ortho

GB

Rematrix Ortho, acellular dermal matrix for soft tissue reconstruction

Product description, indication

Rematrix Ortho is a native, cell-free, non-pyrogenic, resorbable collagen membrane (acellular dermal matrix, ADM) derived from porcine skin intended as a surgical membrane for repair, reinforcement, coverage and reconstruction of soft tissue for supportive use in the repair of tendon and ligament injuries, for example by means of surgical sutures or suture anchors. Examples include

- rotator cuff
- Achilles tendon and other tendons of the ankle joint
- anterior cruciate ligament and other ligaments of the knee joint
- patellar tendon
- tendons of the extensor carpi radialis brevis muscle (tennis elbow)
- biceps tendon
- quadriceps tendon

The acellular collagen matrix is a resilient and biocompatible implant that integrates into the recipient's tissue by means of cellular and microvascular infiltration. While newly regenerated soft tissue is formed in this manner, the implant material is completely resorbed. Rematrix Ortho is ETO sterilised and is intended for single use only. It is supplied in dry form and must be rehydrated before use for at least 5 minutes in physiological saline solution, Ringer's solution or Ringer's lactate solution. It contains no

preservatives and is supplied in a double-barrier sterile packaging. Rematrix Ortho is not artificially cross-linked.

Composition

Sterile, acellular, porcine collagen matrix.

Areas of application

Surgical membrane for repair, reinforcement, coverage and reconstruction of soft tissue for supportive use in the repair of tendon and ligament injuries, for example by means of surgical sutures or suture anchors. Examples:

- rotator cuff
- Achilles tendon and other tendons of the ankle joint
- anterior cruciate ligament and other ligaments of the knee joint
- patellar tendon
- tendons of the extensor carpi radialis brevis muscle (tennis elbow)
- biceps tendon
- quadriceps tendon

Rematrix Ortho is not intended as a replacement of normal body structures or for provision of full mechanical weight-bearing in the repair of tissues such as the rotator cuff, Achilles tendon or other tendons of the ankle joint, the anterior cruciate ligament or other ligaments of the knee joint, the patellar tendon, the biceps or quadriceps tendon or other tendons. The biomechanical load-bearing capacity is provided by the sutures for repairing tendon tears and sutures or suture anchors for attaching the tendon tissue to the bone.

2

Rematrix Ortho

Contraindications

The product is of porcine origin and must not be used in patients with known hypersensitivity to porcine materials.

Potential complications and side effect

The following side effects may result from implantation procedures of acellular matrices: infections, formation of haematomas and seromas, recurrent tissue damage, formation of fistulas, inflammation and formation of adhesions. Since Rematrix Ortho may in some cases be used in surgical wounds whose sterility cannot be guaranteed, antimicrobial agents are often used to help prevent infectious complications. Following such application, wound healing may be slowed down or inhibited. Complications such as delayed wound infection, recurrence of the defect, and the need for repeat surgery must be expected in patients in a critical condition.

Use during pregnancy and breastfeeding

No studies are available on the use of Rematrix Ortho during pregnancy or breastfeeding, nor on its effect on human fertility. The benefits for the mother and the potential risks to the child must therefore be assessed individually by the treating physician before using the product.

Precautions for use

- Do not resterilise Rematrix Ortho.
- Do not use if the packaging has been opened or damaged.
- Do not use either if the seal has been broken or damaged.
- After use, all unused product and packaging materials must be treated and disposed of in accordance with the relevant medical practices and applicable national and local environmental protection laws on the disposal of packaging materials and biological materials.

dance with the relevant medical practices and applicable national and local environmental protection laws on the disposal of packaging materials and biological materials.

- The surgical matrix must not be reused once removed from the packaging and/or after coming into contact with a patient, as this increases the risk of contamination, with a Storage
- Rematrix Ortho is a sterile medical product that must be protected against direct heat and sunlight and stored at room temperature in a clean and dry place.
- Store inside the original packaging.
- The product must not be used after the expiry date.

Sterilisation

This product has been sterilised with ETO and must not be resterilised. Sterility is only guaranteed whilst the packaging is intact.

Application

Rematrix Ortho must only be used by physicians who, through appropriate training, are familiar with tissue regeneration techniques. Since the membrane is of natural origin, a slight corrugation of the collagen structure and a variation of the membrane thickness may occur when dry. This does not affect the quality or functionality of Rematrix Ortho. While ensuring sterility, the product can be cut to the desired size and shape of the defect to be treated either in a dry or rehydrated state using scissors. Appropriate templates can be helpful when determining the required area. The membrane is supplied dry and must be rehydrated prior to implantation by placing it in a sterile dish after opening the packaging

3

Rematrix Ortho

and completely covering it with sterile saline solution, Ringer's solution or Ringer's lactate solution for at least 5 minutes at room temperature. Rematrix Ortho becomes glassy once the rehydration is complete. When implanting the membrane, the greatest possible contact with healthy, well-vascularised tissue must be ensured to promote the ingrowth of cells and the formation of new tissue. An overlap of at least 2-3 cm with healthy tissue is recommended. In case of significant contaminations or infections, in addition to a thorough debridement, measures to control the infection, such as adequate drainage, negative pressure therapy and antimicrobial therapy must be taken prior to implantation. (The presence of a significant bacterial load may compromise the function of the surgical membrane).

Implantation (These recommendations are intended as a general guide only. They take no precedence over institutional protocols or a professional clinical judgement on patient care. Physicians using Rematrix Ortho must be familiar with the surgical procedures and techniques for biological surgical matrices)

1. The site of implantation is to be prepared using standard techniques.
2. While observing aseptic requirements, Rematrix Ortho can be folded and cut to fit the site of implantation while ensuring an overlap.
3. Insert Rematrix Ortho into the implantation site using sterile gloves or sterile forceps.
4. Suture the inserted membrane. Tension and arrangement of the suture depend on the application. A minimum overlap of 2-3 cm with healthy, adjacent tissue provides better results.

5. Complete the standard surgical protocol.
6. Dispose of any unused pieces of the membrane in accordance with the applicable procedures of the facility and applicable environmental regulations.

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