# Rematrix

- GB** Acellular dermal matrix for soft tissue reconstruction  
Instructions for use
- DE** Azelluläre dermale Matrix für Weichgewebe-Rekonstruktion  
Gebrauchsanweisung
- FR** Matrice dermique acellulaire pour la reconstruction de tissus mous  
Notice d'utilisation
- ES** Matriz dérmica acelular para la reconstrucción de tejido blando  
Información de uso
- IT** Matrice dermica acellulare per la ricostruzione di tessuti molli  
Istruzioni per l'uso
- PT** Matriz dérmica acelular para a reconstrução dos tecidos moles  
Informação de uso
- TR** Yumuşak doku rekonstrüksiyonu için hücresiz dermal matriks  
Kullanım talimatları

## Manufacturer

---

medical biomaterial products

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# Rematrix

GB

Acellular dermal matrix for soft tissue reconstruction

---

## Product description

REMATRIX is a cell-free, non-pyrogenic collagen matrix derived from porcine skin (acellular dermal matrix ADM).

REMATRIX is intended for use as a surgical membrane in soft tissue repairs and serves to support, cover or replace tissue. REMATRIX is available in various sizes. REMATRIX comes in double sterile packaging, is supplied dry and does not contain any preservatives.

This surgical membrane is a resilient, biocompatible implant, which is incorporated into the host tissue through cellular and microvascular infiltration.

## Composition

Sterile, acellular, porcine collagen matrix.

## Indications

REMATRIX is intended for implantation in humans. It is suitable for surgical reinforcement and for the protection and/or replacement of soft tissues, promoting tissue regeneration. Today, there are many different indications and possible uses for acellular porcine matrices. A few examples are listed below:

- correction of abdominal and thoracic wall defects,
- treatment of incisional hernias,
- prevention of parastomal hernias.

## Contraindications

The surgical membrane is of porcine origin and cannot be used in patients with known hypersensitivity to porcine materials.

## Warnings

- Do not resterilize.
- Do not use if the packaging has been opened or damaged.
- After use, treat and dispose of all unused parts of the product and all packaging materials in accordance with the recognised medical practices and the applicable national and local environment protection laws on the disposal of packaging materials and biological materials.
- Never reuse the surgical matrix REMATRIX once removed from its packaging and/or after coming into contact with a patient, as this increases the risk of contamination, with subsequent risk of infection.

### **Use during pregnancy and breastfeeding**

There are no researches on the use of REMATRIX in pregnancy and breastfeeding, or its impact on human fertility. Therefore the advantages for the mother and the possible risks for the child must be assessed individually before using REMATRIX

### **Information**

REMATRIX should only be used by physicians with appropriate training in tissue regeneration techniques.

The REMATRIX membrane is of natural origin, thus slight undulating of the collagen structure and variation of the membrane thickness may occur when dry. These effects do not impact the quality and function of REMATRIX.

Ensuring sterility, REMATRIX can be cut to the respective size and shape of the defect to be treated using scissors whilst dry or rehydrated.

The REMATRIX membrane is delivered dry and must be rehydrated before im-

plantation by placing it in a sterile bowl for a minimum of 5 minutes at room temperature after opening the packaging, fully covering it with sterile saline solution, Ringer's solution or Ringer's lactate solution.

During the REMATRIX implantation, ensure the greatest possible contact with well-vascularized tissue to promote cell migration and tissue regeneration. Overlapping with a minimum of 2-3 cm of healthy tissue is then recommended.

In the event of significant contaminations or infections, measures to control the infection must be taken, e.g. suitable drainage, negative pressure therapy and antimicrobial therapy. (The presence of a significant bacterial load could compromise the function of the surgical membrane.)

### **Adverse reactions**

Potential adverse events are those typically associated with acellular matrices and their implantation procedures, including infection, haematoma, seroma

formation, recurring tissue damage, fistula formation, inflammation and adhesion formation.

### Potential complications

Given that REMATRIX may, under certain circumstance, be used in surgical wounds whose sterility is not guaranteed, anti-microbial agents are often used to prevent infectious complications. The use thereof may slow or inhibit wound healing.

In patients in critical condition or with a severely contaminated abdomen, complications such as delayed wound infection, recurrence of the rupture and the need to reoperate must be expected.

### Storage

- REMATRIX is a sterile medical product, which must be stored in a clean and dry location at room temperature, away from direct heat and sunlight.
- Store inside the original packaging.
- Do not use the product after its expiration date.

### Sterilisation

The product is sterilised with ETO and cannot be resterilised. The sterility is only guaranteed whilst the packaging is intact.

### Recommendations for surgical use of REMATRIX

Note: These recommendations merely serve as a general guideline. They do not take precedence over institutional protocols or professional medical judgments on patient care. Physicians must be familiar with the surgical methods and techniques for biological surgical matrices before using REMATRIX.

### Materials required

- Sterile forceps
- Soaking fluid: sterile saline solution, sterile Ringer's solution or sterile Ringer's lactate (room temperature)
- A sterile bowl for each surgical matrix REMATRIX

### Preparation

1. Open the box and remove the double sterile packaging.
2. Open the outer container and remove the inner sterile container using aseptic technique. The inner container is sterile and can be placed directly into the sterile surgical area.
3. Carefully open the inner foil bag and remove the surgical membrane, following the asepsis instructions. Always use sterile gloves or forceps when handling the membrane.
4. Place REMATRIX into a sterile bowl filled with sterile saline solution, sterile Ringer's lactate or sterile Ringer's solution (room temperature) for a minimum of 5 minutes. Make sure that REMATRIX is fully covered by saline solution.
5. Keep REMATRIX in the sterile solution (at room temperature) until ready to use. The acellular matrix can remain in the sterile solution for a maximum of 4 hours.

### Implantation

1. Prepare the implantation site applying the standard techniques.
2. REMATRIX can be folded and cut following asepsis instructions, wherein overlapping with vascularized tissue must be ensured.
3. Insert REMATRIX in the implantation site using sterile gloves or sterile forceps.
4. Suture the membrane.
5. Complete the standard surgical protocol.
6. Dispose of any unused REMATRIX membrane sections following the applicable procedures of the facility and applicable environmental legislation.

Note: Tension and arrangement of the sutures vary by application. A minimum overlap of 2-3 cm with healthy adjacent tissues will provide better results.

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MBP GmbH / Lederstraße 7 / 19306 Neustadt-Glewe