Rematrix Breast

- GB** Acellular dermal matrix for breast reconstruction
Instructions for use
- DE** Azelluläre dermale Matrix für Brustrekonstruktion
Gebrauchsanweisung
- FR** Matrice dermique acellulaire pour la reconstruction mammaire
Notice d'utilisation
- ES** Matriz dérmica acelular para la reconstrucción de mama
Información de uso
- IT** Matrice dermica acellulare per ricostruzione mammaria
Istruzioni per l'uso
- PT** Matriz epidérmica acelular para reconstrução mamária
Informação de uso
- TR** Yumuşak doku rekonstrüksiyonu için hücresiz dermal matriksi
Kullanım talimatları
- PL** Bezkomórkowa macierz skórna do rekonstrukcji tkanek
miękkich piersi
Instrukcja użycia produktu

Manufacturer

mbp

medical biomaterial products

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Rematrix Breast

GB

Acellular dermal matrix for breast reconstruction

Product description

Rematrix Breast is an acellular, non-pyrogenic collagen matrix (acellular dermal matrix ADM), produced from porcine skin. The acellular collagen matrix is a durable and biocompatible implant for soft-tissue reconstruction that is incorporated into the host's tissues by cellular and microvascular infiltration where it is completely absorbed. Rematrix Breast is particularly suitable for one-step breast reconstructions after mastectomy. It is used for the partial or complete cover of silicone grafts that are implanted at the subpectoral or prepectoral (subcutaneous) level. Due to the transformation into peri-implant self tissue, it supports permanent fixation of the implant and achieves good cosmetic outcomes. Rematrix Breast is available in different sizes. Rematrix Breast is supplied dry and contains no preservatives. Rematrix Breast is supplied in a double sterile pack and may only be used once.

Composition

Rematrix Breast is an acellular dermal matrix (ADM) of porcine origin with the main constituent collagen.

Indications

Rematrix Breast is intended for implantation in humans. It is suitable for surgical reinforcement and for the protection and/or replacement of soft tissues, promoting tissue regeneration. Today, there are many different indications and possible uses for acellular porcine matrices. Rematrix Breast is indicated for the one-step breast reconstruction, especially after mastectomy.

Rematrix Breast is indicated for

- subpectoral "one-step" breast reconstruction, especially after mastectomy, and is used for inferior coverage of the lower pole of the implant.
- subcutaneous (prepectoral) "one-step" breast reconstruction and is used for partial or complete coverage of the implant.

Contraindications

Absolute contraindications

The products must not be used for

- patients with known hypersensitivity to porcine material
- patients with tumours infiltrating the skin or chest wall, or that have an

Rematrix Breast

GB

Acellular dermal matrix for breast reconstruction

inflammatory character or a high risk of recurrence.

Relative contraindications

When selecting the patient, it must be noted that the following conditions are associated with a higher risk of complications and require individual consideration of the risk/benefit ratio by the surgeons:

- morbid obesity with an elevated BMI > 40 kg/m²
- inadequately controlled diabetes mellitus
- active smoker
- chronic immunosuppressant therapy
- injury due to previous radiotherapy.

Warnings

- Do not resterilize.
- Do not use if the packaging has been opened or damaged.
- After use, treat and dispose of all unused parts of the product and all packaging materials in accordance with the recognized medical practices and the applicable national and local environment protection laws on the disposal of packaging materials and biological materials.

- Never reuse the surgical matrix once removed from its packaging and/or after coming into contact with a patient, as this increases the risk of contamination, with subsequent risk of infection.

Use during pregnancy and breast-feeding

There are no researches on the use of Rematrix Breast in pregnancy and breast-feeding, or its impact on human fertility. Therefore the advantages for the mother and the possible risks for the child must be assessed individually before using Rematrix Breast.

Information

Rematrix Breast should only be used by doctors who are familiar with techniques of ADM-assisted implant-based breast reconstruction. The Rematrix Breast membrane is of natural origin, thus slight undulating of the collagen structure and variation of the membrane thickness may occur when dry. These effects do not impact the quality and function of Rematrix Breast.

Rematrix Breast

GB

Acellular dermal matrix for breast reconstruction

Ensuring sterility, Rematrix Breast can be cut to fit the desired size and shape in its dry or rehydrated state. The use of suitable phantom can assist to shape the matrix according to the area needed.

The Rematrix Breast membrane is dry at delivery and must be rehydrated before implantation by placing it in a sterile bowl for a minimum of 5 minutes at room temperature after opening the packaging, fully covering it with sterile saline solution, Ringer's solution or Ringer's lactate solution.

During the Rematrix Breast implantation, ensure the greatest possible contact with well-vascularized tissue to promote cell migration and tissue regeneration. Overlapping with a minimum of 2-3 cm of healthy tissue is then recommended.

In the event of significant contaminations or infections, measures to control the infection in accordance with institutional protocols must be taken. (The presence of a significant bacterial load could compromise the function of the surgical membrane.)

Adverse events

The following side effects may occur following implantation procedures with acellular matrices in implant-based breast reconstruction: infections, formation of haematomas and seromas, recurrent tissue damage, fistula formation, inflammation.

Potential complications

Since Rematrix Breast may potentially be implanted in surgical wounds whose sterility is not guaranteed, antimicrobial active substances are often used that are intended to help avoid infectious complications. Wound healing may be slower or inhibited as a result of such use.

The use of acellular matrices to cover breast implants counteracts capsule contractions. However, if Rematrix Breast is used to cover the upper pole of the implant for prepectoral (subcutaneous) breast reconstruction, rippling may occur due to inadequate formation of the subcutis. This can be counteracted by lipomodelling.

Rematrix Breast

GB

Acellular dermal matrix for breast reconstruction

Wound dehiscence, ischaemia of the skin flap and infections must be counteracted by suitable clinical measures; they can lead to loss of the implant.

Storage

- Rematrix Breast is a sterile medical product, which must be stored in a clean and dry location at room temperature, away from direct heat and sunlight.
- Store inside the original packaging.
- Do not use the product after its expiration date.

Sterilization

The product is sterilized with ETO and cannot be resterilized. The sterility is only guaranteed whilst the packaging is intact.

Recommendations for surgical use of Rematrix Breast

Note: These recommendations merely serve as a general guideline. They do not take precedence over institutional protocols or professional medical judgments on patient care. Before using Rematrix Breast, doctors must be familiar with the surgical procedures and techniques

for biological surgical matrices and especially with ADM-assisted implant-based breast reconstruction.

Materials required

- Sterile forceps
- Soaking fluid: sterile saline solution, sterile Ringer's solution or sterile Ringer's lactate (room temperature)
- A sterile bowl for each surgical matrix Rematrix Breast.

Preparation

1. Open the box and remove the double sterile packaging.
2. Open the outer container and remove the inner sterile container using aseptic technique. The inner container is sterile and can be placed directly into the sterile surgical area.
3. Carefully open the inner packaging and remove the surgical membrane, following the asepsis instructions. Always use sterile gloves or forceps when handling the membrane.
4. Place Rematrix Breast into a sterile bowl with sterile saline solution, sterile Ringer's lactate or sterile Ringer's

Rematrix Breast

GB

Acellular dermal matrix for breast reconstruction

solution (room temperature) for a minimum of 5 minutes. Make sure that Rematrix Breast is fully covered by saline solution.

5. Keep Rematrix Breast in the sterile solution (at room temperature) until ready to use. The acellular matrix can remain in the sterile solution for a maximum of 4 hours.

Implantation

1. Prepare the implantation site using standard procedures.
2. Rematrix Breast can be trimmed to suit the intended cover of the breast implant under observance of aseptic technique. If necessary, two membranes can be sutured together.
3. Insert Rematrix Breast into the implantation site with sterile gloves or sterile forceps.
4. Suture the implanted membrane with the disinserted pectoralis muscle and the chest wall (for subpectoral breast reconstruction) or the fascia of the pectoralis muscle (for prepectoral breast reconstruction) around the implanted breast implant.

5. Complete the standard surgical protocol ensuring that there is uninterrupted contact between the membrane and the breast implant and the sutured skin flap.
6. Dispose of the unused parts of the Rematrix Breast membrane under observance of the applicable procedures of the facility and the applicable environmental legislation.

Note: Tension and arrangement of the sutures can vary by application. A minimum overlap of 2-3 cm with healthy adjacent tissues will provide better results.

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MBP GmbH / Lederstraße 7 / 19306 Neustadt-Glewe