medical biomaterial products

# Xenoderm

GB Temporary biological skin replacement  
Instructions for use  
DE Temporärer biologischer Hautersatz  
Gebrauchsanweisung  
FR Remplacement cutané biologique temporaire  
Notice d'utilisation  
ES Sustituto cutáneo biológico temporal  
Información de uso  
IT Trapianto di cute naturale provvisorio  
Istruzioni per l'uso  
PT Substituição biológica temporária da pele  
Informação de uso

REF	UDI	PZN	Size in cm	Content/per Box
0505X	426023090301	2738218	5 x 5	10
1010X	426023090302	2738388	10 x 10	5
1020X	426023090303	2738394	10 x 20	5
1030X	426023090304	2738425	10 x 30	5
1040X	426023090305	-	10 x 40	5
Z0505X	426023090301	-	5 x 5	1
Z1010X	426023090302	-	10 x 10	1
Z1020X	426023090303	-	10 x 20	1
Z1030X	426023090304	-	10 x 30	1
Z1040X	426023090305	-	10 x 40	1

**CE0483**

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# Xenoderm

GB

for Management of skin loss injuries

---

## Definition:

Xenoderm acts as a temporary skin replacement, within the limits of the functions of the patients own skin to:

- Decrease fluid and heat loss
- Enhance the growth of epithelial and granulation tissue
- Protect granulating wounds

## Indications:

- Coverage of full-thickness or partial thickness skin loss due to burns or other trauma
- Coverage of chronic non healing wounds, such as venous, diabetic or pressure ulcers
- Coverage of autograft
- Coverage of donor sites

## Advantages:

- Simple application after reconstitution of Xenoderm
- Good adherence of the biological dressing
- Decrease in pain; exposed nerve endings covered
- No evidence of vascularization within seven days
- No sensitization by the patient to the dressing
- Protection of the underlying tissue against bacterial ingress
- Protection against the drying effects of the environment
- Reapplications are almost painless and without blood loss
- Transparency of the dressing allows inspection of the covered area
- Xenoderm facilitates movement, early ambulation and rehabilitation

## Application:

1. To reconstitute; the inner packet should be opened aseptically and the dressing submerged in a basin containing sterile isotonic saline solution, at room temperature, for at least 2 minutes, or until the dressing becomes translucent.
2. At the time of initial application or reapplication of Xenoderm, the wound area should be thoroughly debrided and gently cleansed, e.g. with an antiseptic; On partial-thickness wounds Xenoderm may be left in place, once adherent, and allowed to slough off as re-epithelialization occurs.
3. Apply Xenoderm firmly to the wound, smoothing to eliminate bubbles and trimming to conform to the wound boundary.
4. As a precaution Xenoderm should not be applied to weight bearing areas of the body (risk of maceration).
5. Xenoderm can be held in place by gauze bandages or any other suitable dressing.
6. Xenoderm should be changed according to the dressing protocol of the unit.
7. Xenoderm dressing changes are almost painless. Always moisten the dressing before removal with sterile isotonic saline solution.
8. The length of time that Xenoderm should be left on the wound is variable and is dictated by the course of wound healing.
9. Infected, heavily secreting wounds may initially require more frequent dressing changes.
10. Xenoderm dressing changes on a daily basis may be required with the above wounds.
11. On all aseptic wounds –including granulating wounds– Xenoderm should be removed within 10 days due to the possibility of excessive adherence.

## Contraindications:

- Full-thickness burn wounds before removal of the eschar
- Badly contaminated wounds which have not been adequately debrided
- Patients with sensitivity to porcine material

## Note:

- Xenoderm is sterilized by  $\gamma$ -irradiation and must not be re-sterilized
- is in a single sterile packet,
  - must be considered non sterile if the packaging is damaged and must be discarded,
  - has a shelf life of five (5) years when stored at room temperature,
  - must not be used after the expiry date,
  - must be protected against moisture,
  - is a medical device and therefore must be kept away from children.

After it has been removed from the packaging and/or has come in contact with a patient, Xenoderm surgical membrane cannot be reused without increased risk of contamination and subsequent infection.

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