Xenoguard

- GB** Pericardium membrane for soft tissue reconstruction
Instructions for use
- DE** Pericardium - Membran für Weichgewebe-Rekonstruktion
Gebrauchsanweisung
- FR** Membrane péricardique pour la reconstruction de tissus mous
Mode d'emploi
- ES** Membrana pericárdica para reconstrucciones de tejidos blandos
Instrucciones de uso
- IT** Membrana pericardica per la ricostruzione di tessuti molli
Istruzioni per l'uso
- PT** Membrana de pericárdio para reconstrução de tecido mole
Informação de uso
- TR** Perikard - Yumuşak doku rekonstrüksiyonu için membran
Kullanım talimatları
- EL** Μεμβράνη από περικάρδιο για την αποκατάσταση ελλειμμάτων
μαλακού ιστού
Οδηγίες χρήσης

Manufacturer

MBP Medical Biomaterial Products GmbH
Lederstraße 7 / 19306 Neustadt-Glewe
info@mbp-gmbh.de / www.mbp-gmbh.de

CE0483

Rev.08/2021-03

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Product description

Xenoguard is a cell-free, non-pyrogenic collagen matrix obtained from porcine pericardium. Xenoguard is intended for use as a surgical membrane in soft tissue repairs and serves to support, cover or replace tissue. Xenoguard is available in various sizes. Xenoguard comes in double sterile packaging, is supplied dry and does not contain any preservatives. This surgical membrane is a resilient, bio-compatible implant, which grows into the host tissue through cellular and microvascular infiltration.

Composition

Sterile, acellular collagen matrix of porcine pericardium.

Indications

Xenoguard is intended for implantation into humans and serves to regenerate soft and connective tissues for

- reconstruction, rebuilding contours and regenerating soft and connective tissues, particularly following loss of tissue,
- supporting tissue in case of loss of stability.

Applications include

- breast reconstruction
- breast and body wall reconstruction
- urological and gynaecological reconstructive surgery of the pelvic floor
- reconstructive ENT surgery

Contraindications

- The surgical membrane is of porcine origin and cannot be used in patients with known hypersensitivity to porcine materials.

Warnings

- Do not resterilise.
- Do not use if the packaging has been opened or damaged. Do not use if the seal is open or damaged.
- After use, treat and dispose of all unused products and packaging materials in accordance with recognised medical practice and the applicable national and local environmental protection laws on the disposal of packaging materials and biological materials.
- Never reuse the surgical matrix Xenoguard once it has been removed from its packaging and/or after it has come

into contact with a patient, as this increases the risk of contamination, with a subsequent risk of infection.

Use during pregnancy and breastfeeding

There are no studies on the use of Xenoguard during pregnancy and breastfeeding or on its impact on human fertility. Therefore the advantages for the mother and the possible risks for the child must be assessed individually before Xenoguard is used.

Application instructions Xenoguard should only be used by physicians who have become familiar with the techniques of tissue regeneration through respective training. As the Xenoguard membrane is of natural origin, slight rippling of the collagen structure and variation of the membrane thickness may occur when it is dry. These effects have no influence on the quality and function of Xenoguard.

While ensuring sterility, Xenoguard can be cut with scissors to the respective size and shape of the defect to be treated in a dry or rehydrated state. The use of suit-

able stencils can assist in determining the area needed.

The Xenoguard membrane is dry on delivery and must be rehydrated before implantation by placing it in a sterile bowl for a minimum of 5 minutes at room temperature after opening the packaging and fully covering it with sterile saline solution, Ringer's solution or Ringer's lactate solution. Xenoguard has a glazed appearance when completely rehydrated.

When inserting the Xenoguard membrane the greatest possible contact with healthy tissue well interspersed with vessels must be ensured in order to promote the ingrowth of cells and the regeneration of tissue. Overlapping with a minimum of 2-3 cm of healthy tissue is recommended.

In the event of significant contamination or infection, in addition to thorough debridement prior to implantation, measures to control the infection must be taken, e.g. suitable drainage, negative pressure therapy and antimicrobial therapy. (The presence of a significant bacter-

ial load could compromise the function of the surgical membrane.)

Adverse reactions

The following adverse reactions can occur in implantation procedures of acellular matrixes: infections, haematoma and seroma formation, recurring tissue damage, fistula formation, inflammation, adhesion formation.

Potential complications

Given that Xenoguard is possibly used in surgical wounds, the sterility of which is not guaranteed, antimicrobial agents are often used to help prevent infectious complications. The use of these agents may slow or inhibit wound healing.

In patients in critical condition or with a severely contaminated abdomen, complications such as delayed wound infection, recurrence of the rupture and the need to reoperate must be expected.

Storage

- Xenoguard is a sterile medical product which must be stored in a clean and

dry location at room temperature away from direct heat and sunlight.

- Store in the original packaging.
- Do not use this product after its expiry date.

Sterilisation

This product was sterilised with ETO and must not be resterilised. The sterility is only guaranteed provided the packaging is intact.

Recommendations for preparing Xenoguard for surgical use

Note: These recommendations merely serve as a general guideline. They do not take precedence over institutional protocols or professional medical judgement on patient care. Physicians must be familiar with the surgical methods and techniques for biological surgical matrixes before using Xenoguard.

Materials required

- Sterile forceps
- Soaking fluid: sterile saline solution, sterile Ringer's solution or sterile Ringer's lactate (room temperature)

Pericardium membrane for soft tissue reconstruction

- A sterile bowl for each piece of the surgical matrix Xenoguard

the sterile solution for a maximum of 4 hours.

Preparation

- Open the box and remove the double sterile packaging.
- Open the outer portion of the sterile packaging and remove the inner portion of the sterile packaging, following the asepsis rules. The inner sterile packaging is sterile and can be introduced directly into the sterile operating area.
- Carefully open the inner foil bag and remove the surgical membrane, following the asepsis rules. Always use sterile gloves or forceps when handling the membrane.
- Place Xenoguard into a sterile bowl filled with sterile saline solution, sterile Ringer's lactate or sterile Ringer's solution (room temperature) for a minimum of 5 minutes. Once fully rehydrated, Xenoguard will have a glazed appearance.
- Keep Xenoguard in the sterile solution (at room temperature) until ready to use. The acellular matrix can remain in

Implantation

- Prepare the implantation site using the standard techniques.
- Xenoguard can be folded and cut to size following the asepsis rules, whereby overlapping must be ensured.
- Insert Xenoguard into the implantation site using sterile gloves or sterile forceps.
- Suture the inserted membrane.
- Complete the standard surgical protocol.
- Dispose of any unused Xenoguard membrane parts following the applicable procedures of the facility and the applicable environmental legislation.

Note: The tension and arrangement of the sutures depend on the application. When repairing gaps, surgical experience with soft tissue implants has shown that suturing the surgical matrix Xenoguard under tension while overlapping the healthy adjacent tissue by a minimum of 2-3 cm will provide better results.

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